

Cold Chain Management Solution

Simple Installation with Drastic Improvements



Cold Chain Management Solution: Simple Installation with Drastic Improvements

Designed for chemists, pharmacists and anyone who needs to manage refrigerated products for the health sector.

A simple to install, wireless IoT solution ensuring the safety, efficacy and compliance of stored vaccines, pharmaceuticals and other temperature sensitive products.



NONASH University



Monitoring Covid-19 vaccines right now at **Barwon Health's University Hospital** and biospecimens at **Monash BioMedicine Discovery Institute (BDI), St Vincent's Hospital** and **NSW HealthShare** facilities.



How it Works

The CT Cold Chain Management Solution is a cloud-based software platform that pairs with our latest IoT wireless smart tag to comprehensively monitor product temperatures, ensuring product safety, as well as providing operational and equipment insights.

CT's IoT smart tags are calibrated according to NATA standards and can monitor extreme temperatures to the requisite accuracy for highly sensitive materials such as vaccines.

With our cloud-based software platform and mobile application, we deliver the ability to convert your paper-based processes to digital forms, providing you with a complete end to end compliance solution.



COLD CHAIN MANAGEMENT SOLUTION

Comprehensive Dashboards and Mobile application to visualise your product compliance from anywhere.


Real Time Dashboard					Real Time	Dashboard				
Shipment Dashboard	Export to Excel	Export to PDF							Last Refreshed at 3	1/05/2021 11:38
Locations	Locations	Description	Туре	Value	Status	Signal	Battery	Alert	Graph	Notes
	Ŧ	T		T						
Profiles										
Users	112041 - Robot 1 Room 10600390			24.6 °C	(It-	atil	—	8	1~~.	t.
Reports	114361 - Robot 2 Room 10600560		8	26.1 °C	(It-	at		8	Im.	t.
Gateway	111000 - Pharmacy Freezer 1	BELL - 03	8	-20.7 °C	((i-	all.		8	1~	C.
Sites Add Devices	111016 - Pharmacy Fridge 1	BELL - 03		3.8 °C	((i·	att		8	1~	C.
Add Gateways	ALCC Freezer	ALCC - 03	8	-20.5 °C	((i·	all	—	8	i~	C.
Corrective Actions	L6 SJOG Ophthalmology Drug Fridge	SJOG LVL 6	8	5.5 °C	((r·	aff			1	C.
Alert Reasons Manage Shipment	111013 - SJOG Pharmacy Drug Fridge	SJOG Lvl 1	8	5.4 °C	((t·	all.			I~	Ľ⊕
Assign Form	111038 - McKellar Pharmacy	MCKIRC - 01	8	4.3 °C	(((-	all	—		1~	C.
Site Forms	109325 - Satellite Drug Fridge	CC - 06	8	5.1 °C	((i·	all			1~~	C.
Forms Report	TBA - Pathology							~		
Staff List	н 4 1 2 🕨	н								1 - 20 of 3
Manage Form Templates										

SITE OVERVIEW

COLD CHAIN MANAGEMENT SOLUTION

Real Time Dashboard			Location Na	me: 109325 - Satellite	Drug Fridge	14D 30D	×		
Shipment Dashboard	Export to Excel		on		30 10			Last Refreshed	at 31/05/2021 12:05 💭
Shipment Dashboard	Locations			Temperature		S	± _	Graph	Notes
Locations									
Profiles		8.00	Temperature Lower To	lerance Opper Tolerance	Amber Lower Toleranc	Amber Upper Toleranc	<u>.</u>		
Users	112041 - Robot Room 1060039	e.00		and the second second second	htter and the second	****	WW	<u>I~~</u>	C
Reports	114361 - Robot Room 1060056	ad ∎ 4.00						1~~	Ľ⊕
Gateway	111000 - Pharma Freezer 1	2.00 28/05/21 12:10:34	29/05/21 11:00:00	30/05/21 00:53:20	30/05/21 14:46:40	31/05/21 08:	50:34	I~.	C
Sites	111016 - Pharma	28/05/21 12	29/00/**	Sample Time	30/09/21	31/05/21			
Add Devices	Fridge 1						CLOSE	I~	Ľ⊕
Add Gateways	ALCC Freezer	8	-20.0-0	.Ċ.	IIIa	<u></u>	8	1~	C⊕
Corrective Actions	L6 SJOG Ophthalmology Drug Fridge	8	5.1 °C	(((-	att			1~~.	C⊕
Alert Reasons	111013 - SJOG		5.4 °C	()	100		R	ter.	C.
Manage Shipment	Pharmacy Drug Fridge	, •	5.4 °C	1	100	<u></u>	X	~	L⊕
Assign Form	111038 - McKellar Pharmacy	8	4.3 °C	((t·	and the second			<u>Inc</u>	Ľ⊕
Site Forms	109325 - Satellite Drug Fridge	8	5.1 °C	(((:	all			1~~	t⊕
Forms Report	TBA - Pathology specimen storage		49 °C	6	all		8	tw	۲۳.
Staff List	₩ 4 1 2 ►	н							1 - 20 of 33
Manage Form Templates									

DRILL DOWN INTO INDIVIDUAL ASSETS



SMART ANALYTICS

Key Features

Accurate costeffective monitoring at your fingertips with a simple monthly subscription. Receive near real-time alarms to identify potential issues and take proactive action. Generate reports for compliance purposes.



Easy to Install

CT's IoT Smart Tags are wireless battery-powered devices, designed to be easy to install. Simply turn it on and place it within your existing refrigeration assets to provide continuous automated temperature monitoring.



Compliance Alerts

Benefit from temperature threshold alerts sent to you via SMS and email, and record any corrective actions. Setup manual temperature monitoring reminders to meet Strive for 5 Compliance manual monitoring requirements.



CT Compliance Management

Migrate from paper-based forms to digital. Our solution allows you to create unlimited digital forms to help you centralise and manage your vaccine compliances and other Pharmacy task management requirements.

A.				
d_	0			
-	ь			

Smart Analytics

CT's advanced analytics and alerts provide insights to empower your organisation to proactively detect refrigeration asset failures, minimise maintenance costs and reduce product wastage.



Visualisation & Reporting

Dynamic visualisation and reporting platform for desktop and mobile, allowing you to easily share information with key stakeholders, and generate data for any defined time period.



Secure & Encrypted

Benefit from CT's proprietary encryption algorithm and protocols to ensure your data is safe and secure. We take data sovereignty seriously, and only host data within Australian datacentres.

Talk to our Cold Chain Solution Specialist.

\$ 1800 100 227



Level 7, 420 Collins Street Melbourne Victoria 3000 Australia

Tel 1800 100 227 Email sales@ct1limited.com

constellationtechnologies.com.au